



Urine Monovette® User Guide




Remove stopper and keep for later use!


Attach plastic straw.



Insert the plastic straw into the container and fill the Urine Monovette® up to the bottom of the label.



To empty the plastic tube, hold the Urine Monovette® in an upright position and pull the plunger backwards to the bottom of the tube.



Remove the plastic straw break off the plunger and throw away.
Replace the stopper.

Technical modifications reserved



SARSTEDT

SARSTEDT Ltd. · 68 Boston Rd, Leicester LE4 1AW · Tel.: 0116 235 9023 · Fax: 0116 236 6099

info@sarstedt.co.uk · www.sarstedt.com